

# NDC80 antibody



Catalog Number: 113052

## Product name

NDC80 antibody

IP: 1:200-1:2000

IHC: 1:20-1:200


## Specificity

Human, Mouse, Rat; other species not tested.

## Validations

## Antibody description

NDC80 Rabbit Polyclonal antibody. Positive IHC detected in human brain tissue. Positive WB detected in A2780 cells, human spleen tissue. Positive IP detected in A2780 cells. Observed molecular weight by Western-blot: 76 kDa




## Preparation

This antibody was obtained by immunization of NDC80 recombinant protein (Accession Number: NM\_006101). Purification method: Antigen affinity purified.

A2780 cells were subjected to SDS PAGE followed by western blot with Catalog No:113052(NDC80 antibody) at dilution of 1:500

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT


## Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human brain using Catalog No:113052(NDC80 antibody) at dilution of 1:100 (under 10x lens)

## Ig Type

Rabbit IgG



## Applications

ELISA, WB, IHC, IP

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

IP Result of anti-NDC80 (IP:Catalog No:113052, 4ug; Detection:Catalog No:113052 1:500) with

# NDC80 antibody



Catalog Number: 113052

---

A2780 cells lysate 960ug.