

NDRG2 antibody



Catalog Number: 113057

Product name

NDRG2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NDRG2 Rabbit Polyclonal antibody. Positive WB detected in rat brain tissue, mouse brain tissue. Positive IP detected in rat brain tissue. Observed molecular weight by Western-blot: 41-43 kDa

Preparation

This antibody was obtained by immunization of NDRG2 recombinant protein (Accession Number: XM_011537001). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:500-1:5000

Validations



rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113057(NDRG2 antibody) at dilution of 1:1000



IP Result of anti-NDRG2 (IP:Catalog No:113057, 4ug; Detection:Catalog No:113057 1:1000) with rat brain tissue lysate 4000ug.