

# NDRG2 antibody

Catalog Number: 113057

## Product name

NDRG2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

NDRG2 Rabbit Polyclonal antibody. Positive WB detected in rat brain tissue, mouse brain tissue. Positive IP detected in rat brain tissue. Observed molecular weight by Western-blot: 41-43 kDa

## Preparation

This antibody was obtained by immunization of NDRG2 recombinant protein (Accession Number: XM\_011537001). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:500-1:5000

## Validations



rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113057(NDRG2 antibody) at dilution of 1:1000



IP Result of anti-NDRG2 (IP:Catalog No:113057, 4ug; Detection:Catalog No:113057 1:1000) with rat brain tissue lysate 4000ug.