

NDUFA8 antibody



Catalog Number: 113064

Product name

NDUFA8 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

NDUFA8 Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, mouse heart tissue, mouse liver tissue. Observed molecular weight by Western-blot: 19kd

Preparation

This antibody was obtained by immunization of NDUFA8 recombinant protein (Accession Number: NM_014222). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:113064(NDUFA8 antibody) at dilution of 1:500