

NDUFB3 antibody



Catalog Number: 113073

Product name

NDUFB3 antibody

Immunogen

[Human NDUFB3 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NDUFB3 Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, mouse heart tissue, rat heart tissue. Observed molecular weight by Western-blot: 11 kDa

Preparation

This antibody was obtained by immunization of NDUFB3 recombinant protein (Accession Number: NM_001257102). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:113073(NDUFB3 antibody) at dilution of 1:500