

# NDUFB5 antibody



Catalog Number: 113074

## Product name

NDUFB5 antibody

## Immunogen

[Human NDUFB5 Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

NDUFB5 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, mouse spleen tissue. Observed molecular weight by Western-blot: 16 kDa

## Preparation

This antibody was obtained by immunization of NDUFB5 recombinant protein (Accession Number: NM\_002492). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:113074(NDUFB5 Antibody) at dilution of 1:600