

NDUFB9 antibody



Catalog Number: 113078

Product name

NDUFB9 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

NDUFB9 Rabbit Polyclonal antibody. Positive WB detected in mouse skin tissue, A375 cells, human liver tissue, human pancreas tissue, mouse pancreas tissue. Observed molecular weight by Western-blot: 22kd

Preparation

This antibody was obtained by immunization of NDUFB9 recombinant protein (Accession Number: NM_005005). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse skin tissue were subjected to SDS PAGE followed by western blot with Catalog No:113078(NDUFB9 antibody) at dilution of 1:500