

NDUFS7 antibody



Catalog Number: 113085

Product name

NDUFS7 antibody

IHC: 1:20-1:200

IF: 1:10-1:100

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NDUFS7 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, K-562 cells. Positive IHC detected in human liver tissue. Positive IF detected in Hela cells. Observed molecular weight by Western-blot: 20kd

Preparation

This antibody was obtained by immunization of NDUFS7 recombinant protein (Accession Number: NM_024407). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:1000-1:10000


Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113085(NDUFS7 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human liver using Catalog No:113085(NDUFS7 antibody) at dilution of 1:100 (under 10x lens)




Immunohistochemical of paraffin-embedded human liver using Catalog No:113085(NDUFS7

NDUFS7 antibody



Catalog Number: 113085

antibody) at dilution of 1:100 (under 40x lens)



Immunofluorescent analysis of Hela cells, using NDUFS7 antibody Catalog No:113085 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).