

# NEDD8 antibody

Catalog Number: 113098

## Product name

NEDD8 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

NEDD8 Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue. Positive WB detected in mouse brain tissue, human heart tissue. Observed molecular weight by Western-blot: 6-10 kDa

## Preparation

This antibody was obtained by immunization of NEDD8 recombinant protein (Accession Number: NM\_006156). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC

## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113098(NEDD8 antibody) at dilution of 1:800



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using Catalog No:113098(NEDD8 Antibody) at dilution of 1:50 (under 40x lens)