

NEIL3 antibody



Catalog Number: 113100

Product name

NEIL3 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

NEIL3 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, A431 cells, human skin tissue, mouse skin tissue. Observed molecular weight by Western-blot: 62-68 kDa

Preparation

This antibody was obtained by immunization of NEIL3 recombinant protein (Accession Number: NM_018248). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:113100(NEIL3 antibody) at dilution of 1:500