

NEK11 antibody

Catalog Number: 113101

Product name

NEK11 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NEK11 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HepG2 cells. Observed molecular weight by Western-blot: 54kd

Preparation

This antibody was obtained by immunization of NEK11 recombinant protein (Accession Number: NM_001321223). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113101(NEK11 antibody) at dilution of 1:400