

NEK9 antibody

Catalog Number: 113107

Product name

NEK9 antibody

Specificity

Human; other species not tested.

Antibody description

NEK9 Rabbit Polyclonal antibody. Positive IF detected in COS-7 cells.

Preparation

This antibody was obtained by immunization of NEK9 recombinant protein (Accession Number: BC009336). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IF

Dilutions

Recommended Dilution:

IF: 1:20-1:200

Validations



Immunofluorescent analysis of (10% Formaldehyde) fixed COS-7 cells using Catalog No:113107(NEK9 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)