

Netrin G1 antibody



Catalog Number: 113118

Product name

Netrin G1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Netrin G1 Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue, A549 cells, human brain tissue, human testis tissue, mouse brain tissue. Observed molecular weight by Western-blot: 60 kDa

Preparation

This antibody was obtained by immunization of Netrin G1 recombinant protein (Accession Number: NM_001372168). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:113118(NTNG1 antibody) at dilution of 1:500