

# NEURL3 antibody

Catalog Number: 113123

## Product name

NEURL3 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

NEURL3 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 32 kDa

## Preparation

This antibody was obtained by immunization of NEURL3 recombinant protein (Accession Number: BC012317). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IP

## Dilutions

Recommended Dilution:

IP: 1:200-1:1000

## Validations



IP Result of anti-NEURL3 (IP:Catalog No:113123, 3ug; Detection:Catalog No:113123 1:300) with mouse brain tissue lysate 5000ug.