

NEUROD4 antibody



Catalog Number: 113126

Product name

NEUROD4 antibody

Polyclonal

Ig Type

Specificity

Human, Mouse, Rat; other species not tested.

Rabbit IgG

Applications

Antibody description

NEUROD4 Rabbit Polyclonal antibody. Positive WB detected in Y79 cells, Neuro-2a cells. Observed molecular weight by Western-blot: 50-55 kDa

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of NEUROD4 recombinant protein (Accession Number: BC040961). Purification method: Antigen affinity purified.

Formulation


PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Validations



Y79 cells were subjected to SDS PAGE followed by western blot with Catalog No:113126(NEUROD4 Antibody) at dilution of 1:1000