

# neuromedin B antibody



Catalog Number: 113130

---

## Product name

neuromedin B antibody

ELISA, IHC

## Dilutions

Recommended Dilution:

IHC: 1:10-1:100

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

neuromedin B Rabbit Polyclonal antibody. Positive IHC detected in human brain tissue, human breast cancer tissue, human skin tissue, human stomach tissue.

## Preparation

This antibody was obtained by immunization of neuromedin B recombinant protein (Accession Number: NM\_205858). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality


Polyclonal

## Ig Type


Rabbit IgG

## Applications

## Validations



Immunohistochemical of paraffin-embedded human brain using Catalog No:113130(NMB antibody) at dilution of 1:25 (under 10x lens)



Immunohistochemical of paraffin-embedded human brain using Catalog No:113130(NMB antibody) at dilution of 1:25 (under 40x lens)