

neurturin-Specific antibody



Catalog Number: 113135

Product name

neurturin-Specific antibody

Specificity

Human; other species not tested.

Antibody description

neurturin-Specific Rabbit Polyclonal antibody.
Positive IHC detected in human colon tissue.
Positive WB detected in human liver tissue, human heart tissue. Observed molecular weight by Western-blot: 22 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_004558).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:113135(NRTN-Specific antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human colon using Catalog No:113135(NRTN-Specific antibody) at dilution of 1:50 (under 10x lens)