

NFE2L3 antibody



Catalog Number: 113147

Product name

NFE2L3 antibody

Immunogen

Specificity

human; other species not tested.

Antibody description

NFE2L3 Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue. Observed molecular weight by Western-blot: 76kd

Preparation

This antibody was obtained by immunization of NFE2L3 recombinant protein (Accession Number: NM_004289). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:113147(NFE2L3 antibody) at dilution of 1:500