

# NFIB antibody



Catalog Number: 113151

---

## Product name

NFIB antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

NFIB Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, HepG2 cells, L02 cells. Observed molecular weight by Western-blot: 51-55kd

## Preparation

This antibody was obtained by immunization of NFIB recombinant protein (Accession Number: NM\_005596). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:113151(NFIB antibody) at dilution of 1:500