

# NFKBIL1 antibody



Catalog Number: 113160

## Product name

NFKBIL1 antibody

Polyclonal

## Ig Type

Rabbit IgG

## Specificity

Human; other species not tested.

## Applications

ELISA, WB

## Antibody description

NFKBIL1 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, MDA-MB-453s cells. Observed molecular weight by Western-blot: 40-43 kDa

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000

## Preparation

This antibody was obtained by immunization of NFKBIL1 recombinant protein (Accession Number: NM\_005007). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:113160(NFKBIL1 Antibody) at dilution of 1:1000

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality