

NFU1 antibody

Catalog Number: 113165

Product name

NFU1 antibody

IP: 1:500-1:5000

IF: 1:10-1:100

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NFU1 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive IP detected in mouse brain tissue. Positive WB detected in mouse cerebellum tissue, mouse brain tissue, mouse testis tissue. Observed molecular weight by Western-blot: 23-28 kDa

Preparation

This antibody was obtained by immunization of NFU1 recombinant protein (Accession Number: NM_001002755). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IF, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



mouse cerebellum tissue were subjected to SDS PAGE followed by western blot with Catalog No:113165(NFU1 Antibody) at dilution of 1:1000



IP Result of anti-NFU1 (IP:Catalog No:113165, 4ug; Detection:Catalog No:113165 1:1000) with mouse brain tissue lysate 4000ug.



Immunofluorescent analysis of HepG2 cells using Catalog No:113165(NFU1 Antibody) at dilution of

NFU1 antibody



Catalog Number: 113165

1:25 and Alexa Fluor 488-conjugated AffiniPure

Goat Anti-Rabbit IgG(H+L)