

NFYA antibody

Catalog Number: 113167

Product name

NFYA antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NFYA Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in HEK-293 cells, HeLa cells, HT-1080 cells, Jurkat cells, K-562 cells. Observed molecular weight by Western-blot: 40kd, 46kd

Preparation

This antibody was obtained by immunization of NFYA recombinant protein (Accession Number: NM_021705). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:113167(NFYA antibody) at dilution of 1:200



Immunofluorescent analysis of HepG2 cells, using NFYA antibody Catalog No:113167 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).