

NFYC antibody



Catalog Number: 113168

Product name

NFYC antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NFYC Rabbit Polyclonal antibody. Positive WB detected in K-562 cells. Observed molecular weight by Western-blot: 48 kDa

Preparation

This antibody was obtained by immunization of NFYC recombinant protein (Accession Number: XM_017001368). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:113168(NFYC antibody) at dilution of 1:1000