

# NGEF antibody



Catalog Number: 113172

## Product name

NGEF antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

NGEF Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue, mouse testis tissue. Observed molecular weight by Western-blot: 75 kDa

## Preparation

This antibody was obtained by immunization of NGEF recombinant protein (Accession Number: XM\_011510923). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113172(NGEF antibody) at dilution of 1:500



IP Result of anti-NGEF (IP:Catalog No:113172, 4ug; Detection:Catalog No:113172 1:500) with mouse brain tissue lysate 4000ug.