

NGFRAP1 antibody



Catalog Number: 113173

Product name

NGFRAP1 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

NGFRAP1 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue. Observed molecular weight by Western-blot: 28kd

Preparation

This antibody was obtained by immunization of NGFRAP1 recombinant protein (Accession Number: NM_014380). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:113173(NGFRAP1 antibody) at dilution of 1:400