

# NKX3-1 antibody



Catalog Number: 113205

## Product name

NKX3-1 antibody

Polyclonal

## Ig Type

Rabbit IgG

## Specificity

human, mouse, rat; other species not tested.

## Applications

ELISA, WB

## Antibody description

NKX3-1 Rabbit Polyclonal antibody. Positive WB detected in human testis tissue. Observed molecular weight by Western-blot: 36kd

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000

## Preparation

This antibody was obtained by immunization of NKX3-1 recombinant protein (Accession Number: NM\_006167). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:113205(NKX3-1 antibody) at dilution of 1:800

## Clonality