

NPEPL1 antibody



Catalog Number: 113211

Product name

NPEPL1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NPEPL1 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in HepG2 cells, HeLa cells, MCF7 cells. Observed molecular weight by Western-blot: 56 kDa

Preparation

This antibody was obtained by immunization of NPEPL1 recombinant protein (Accession Number: NM_024663). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:113211(NPEPL1 antibody) at dilution of 1:50



Immunofluorescent analysis of (10% Formaldehyde) fixed HepG2 cells using Catalog No:113211(NPEPL1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)