

NPHS2 antibody



Catalog Number: 113216

Product name

NPHS2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NPHS2 Rabbit Polyclonal antibody. Positive IF detected in zebrafish embryos. Positive WB detected in mouse colon tissue, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 34kd

Preparation

This antibody was obtained by immunization of NPHS2 recombinant protein (Accession Number: NM_001297575). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:20-1:200

Validations



The podocin antibody Catalog No:113216 recognized the zebrafish pronephric glomeruli in zebrafish embryos(72 hour-old fish). IF result from



Dr. Weibin Zhou.

mouse colon tissue were subjected to SDS PAGE followed by western blot with Catalog No:113216(NPHS2 antibody) at dilution of 1:400