

# NPT1 antibody



Catalog Number: 113223

---

## Product name

NPT1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

NPT1 Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue, mouse kidney tissue. Observed molecular weight by Western-blot: 50 kDa, 60-64 kDa

## Preparation

This antibody was obtained by immunization of NPT1 recombinant protein (Accession Number: XM\_017011201). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:113223(NPT1 antibody) at dilution of 1:300