

# NR1H2 antibody



Catalog Number: 113231

---

## Product name

NR1H2 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

NR1H2 Rabbit Polyclonal antibody. Positive WB detected in human liver tissue, mouse liver tissue. Observed molecular weight by Western-blot: 50 kDa, 40 kDa

## Preparation

This antibody was obtained by immunization of NR1H2 recombinant protein (Accession Number: NM\_007121). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:113231(NR1H2 antibody) at dilution of 1:500