

NR1H3 antibody



Catalog Number: 113232

Product name

NR1H3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NR1H3 Rabbit Polyclonal antibody. Positive WB detected in Sp2/0 cells, rat colon tissue. Observed molecular weight by Western-blot: 44-48 kDa

Preparation

This antibody was obtained by immunization of NR1H3 recombinant protein (Accession Number: XM_047426198). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Sp2/0 cells were subjected to SDS PAGE followed by western blot with Catalog No:113232(NR1H3 antibody) at dilution of 1:800