

NLGN4Y antibody



Catalog Number: 113239

Product name

NLGN4Y antibody

Specificity

Human; other species not tested.

Antibody description

NLGN4Y Rabbit Polyclonal antibody. Positive WB detected in SGC-7901 cells. Observed molecular weight by Western-blot: 72 kDa

Preparation

This antibody was obtained by immunization of NLGN4Y recombinant protein (Accession Number: BC032567). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



SGC-7901 cells were subjected to SDS PAGE followed by western blot with Catalog No:113239(NLGN4Y antibody) at dilution of 1:600