

# NLRP1 antibody

Catalog Number: 113241

## Product name

NLRP1 antibody

## Specificity

Human; other species not tested.

## Antibody description

NLRP1 Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer tissue. Positive WB detected in Transfected HEK-293 cells. Observed molecular weight by Western-blot: 160-166kd

## Preparation

This antibody was obtained by immunization of NLRP1 recombinant protein (Accession Number: BC016965). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:113241(NLRP1 antibody) at dilution of 1:50 (under 10x lens)



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:113241(NLRP1 Antibody) at dilution of 1:1000