

# NMDAR2A antibody



Catalog Number: 113245

## Product name

NMDAR2A antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

NMDAR2A Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Positive IF detected in Hela cells. Observed molecular weight by Western-blot: 150-160 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_000833).  
Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113245 (GRIN2A antibody) at dilution of 1:500



Immunofluorescent analysis of Hela cells, using GRIN2A antibody Catalog No:113245 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).