

# NMNAT1 antibody



Catalog Number: 113252

## Product name

NMNAT1 antibody

## Specificity

human; other species not tested.

## Antibody description

NMNAT1 Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer. Positive WB detected in human testis tissue, human brain tissue, human heart tissue. Observed molecular weight by Western-blot: 32 kDa

## Preparation

This antibody was obtained by immunization of NMNAT1 recombinant protein (Accession Number: XM\_047428076). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:113252(NMNAT1 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:113252(NMNAT1 antibody) at dilution of 1:50 (under 10x lens)