

# NMNAT3 antibody



Catalog Number: 113253

## Product name

Polyclonal

NMNAT3 antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

## Applications

## Antibody description

ELISA, WB

NMNAT3 Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 24-28 kDa

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Preparation

This antibody was obtained by immunization of NMNAT3 recombinant protein (Accession Number: XM\_017006338). Purification method: Antigen affinity purified.

## Formulation


PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

## Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:113253(NMNAT3 Antibody) at dilution of 1:300