

Product name

nNOS antibody

Specificity

Human; other species not tested.

Antibody description

nNOS Rabbit Polyclonal antibody. Positive IHC detected in human kidney tissue.

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_000620). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, IHC


Dilutions

Recommended Dilution:

IHC: 1:20-1:200

Validations

Immunohistochemical of paraffin-embedded human kidney using Catalog No:113259(nNOS antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human kidney using Catalog No:113259(nNOS antibody) at dilution of 1:100 (under 40x lens)