

# NOB1 antibody

Catalog Number: 113262

## Product name

NOB1 antibody

IP: 1:200-1:1000

IHC: 1:10-1:100

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

NOB1 Rabbit Polyclonal antibody. Positive IHC detected in human lymphoma tissue, human colon cancer tissue. Positive IP detected in HepG2 cells. Positive WB detected in SGC-7901 cells. Observed molecular weight by Western-blot: 47-50 kDa

## Preparation

This antibody was obtained by immunization of NOB1 recombinant protein (Accession Number: BC000050). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IHC, IP

## Dilutions


Recommended Dilution:

WB: 1:200-1:1000


## Validations



SGC-7901 cells were subjected to SDS PAGE followed by western blot with Catalog No:113262(NOB1 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:113262(NOB1 antibody) at dilution of 1:25 (under 10x lens)



IP Result of anti-NOB1 (IP:Catalog No:113262, 4ug; Detection:Catalog No:113262 1:300) with

# NOB1 antibody



Catalog Number: 113262

---

HepG2 cells lysate 2400ug.