

NOB1 antibody



Catalog Number: 113262

Product name

NOB1 antibody

IP: 1:200-1:1000

IHC: 1:10-1:100

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NOB1 Rabbit Polyclonal antibody. Positive IHC detected in human lymphoma tissue, human colon cancer tissue. Positive IP detected in HepG2 cells. Positive WB detected in SGC-7901 cells. Observed molecular weight by Western-blot: 47-50 kDa

Preparation

This antibody was obtained by immunization of NOB1 recombinant protein (Accession Number: BC000050). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:1000


Validations



SGC-7901 cells were subjected to SDS PAGE followed by western blot with Catalog No:113262(NOB1 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:113262(NOB1 antibody) at dilution of 1:25 (under 10x lens)



IP Result of anti-NOB1 (IP:Catalog No:113262, 4ug; Detection:Catalog No:113262 1:300) with

NOB1 antibody



Catalog Number: 113262

HepG2 cells lysate 2400ug.