

# NR2C2 antibody



Catalog Number: 113265

## Product name

NR2C2 antibody

## Specificity

Human; other species not tested.

## Antibody description

NR2C2 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, Jurkat cells, PC-3 cells. Observed molecular weight by Western-blot: 70 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_003298). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113265(NR2C2 Antibody) at dilution of 1:1500