

NR2C2 antibody



Catalog Number: 113265

Product name

NR2C2 antibody

Specificity

Human; other species not tested.

Antibody description

NR2C2 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, Jurkat cells, PC-3 cells. Observed molecular weight by Western-blot: 70 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_003298). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113265(NR2C2 Antibody) at dilution of 1:1500