

# NR5A2 antibody



Catalog Number: 113270

## Product name

NR5A2 antibody

Polyclonal

## Ig Type

Rabbit IgG

## Specificity

Human, Mouse, Rat; other species not tested.

## Applications

ELISA, WB

## Antibody description

NR5A2 Rabbit Polyclonal antibody. Positive WB detected in mouse ovary tissue, rat testis tissue. Observed molecular weight by Western-blot: 61-67 kDa

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000

## Preparation

This antibody was obtained by immunization of NR5A2 recombinant protein (Accession Number: NM\_003822). Purification method: Antigen Affinity purified.

## Validations



## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

mouse ovary tissue were subjected to SDS PAGE followed by western blot with Catalog No:113270(NR5A2 Antibody) at dilution of 1:1000

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality