

# NRBP2 antibody

Catalog Number: 113273

## Product name

NRBP2 antibody

## Immunogen

[Human NRBP2 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

NRBP2 Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, HeLa cells, human placenta tissue, human testis tissue, mouse brain tissue. Positive IP detected in mouse kidney tissue. Observed molecular weight by Western-blot: 33 kDa, 58 kDa

## Preparation

This antibody was obtained by immunization of NRBP2 recombinant protein (Accession Number: BC113873). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:113273(NRBP2 antibody) at dilution of 1:500



IP Result of anti-NRBP2 (IP:Catalog No:113273, 4ug; Detection:Catalog No:113273 1:500) with mouse kidney tissue lysate 4000ug.