

NRD1 antibody

Catalog Number: 113275

Product name

NRD1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NRD1 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells. Positive IHC detected in human testis tissue. Observed molecular weight by Western-blot: 135-140kd

Preparation

This antibody was obtained by immunization of NRD1 recombinant protein (Accession Number: NM_001101662). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:113275(NRD1 antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human testis using Catalog No:113275(NRD1 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human testis using Catalog No:113275(NRD1 antibody) at dilution of 1:50 (under 40x lens)