

NRXN1 antibody



Catalog Number: 113283

Product name

NRXN1 antibody

Specificity

Human; other species not tested.

Antibody description

NRXN1 Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells. Observed molecular weight by Western-blot: 160-170 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_004801). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:113283(NRXN1 antibody) at dilution of 1:500