

# NOD2 antibody



Catalog Number: 113290

---

## Product name

NOD2 antibody

## Specificity

human, mouse; other species not tested.

## Antibody description

NOD2 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, COLO 320 cells, U-937 cells, RAW264.7 cells. Observed molecular weight by Western-blot: 100-110kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_022162). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:113290(NOD2 antibody) at dilution of 1:800