

# NOL8 antibody

Catalog Number: 113297

## Product name

NOL8 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

NOL8 Rabbit Polyclonal antibody. Positive IF detected in HepG2. Positive WB detected in SGC-7901 cells. Observed molecular weight by Western-blot: 150kd

## Preparation

This antibody was obtained by immunization of NOL8 recombinant protein (Accession Number: BC013788). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

## Validations



SGC-7901 cells were subjected to SDS PAGE followed by western blot with Catalog No:113297(NOL8 antibody) at dilution of 1:500



Immunofluorescent analysis of HepG2 cells, using NOL8 antibody Catalog No:113297 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).