

NOP56 antibody



Catalog Number: 113304

Product name

NOP56 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

NOP56 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, A549 cells. Observed molecular weight by Western-blot: 58-66kd

Preparation

This antibody was obtained by immunization of NOP56 recombinant protein (Accession Number: NM_006392). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113304(NOP56 antibody) at dilution of 1:500