

# Notch1 antibody

Catalog Number: 113307

## Product name

Notch1 antibody

## Specificity

Human; other species not tested.

## Antibody description

Notch1 Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 110-120 kDa

## Preparation

This antibody was obtained by immunization of Notch1 recombinant protein (Accession Number: NM\_008714). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113307(Notch1 antibody) at dilution of 1:500



IP Result of anti-Notch1 (IP:Catalog No:113307, 5ug; Detection:Catalog No:113307 1:200) with HeLa cells lysate 2000ug.