

# NOX2 antibody



Catalog Number: 113312

## Product name

NOX2 antibody

## Immunogen

[Human CYBB Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

NOX2 Rabbit Polyclonal antibody. Positive IHC detected in human skin tissue. Positive WB detected in mouse brain tissue, human serum tissue, mouse heart tissue, Raji cells. Observed molecular weight by Western-blot: 55 kDa

## Preparation

This antibody was obtained by immunization of NOX2 recombinant protein (Accession Number: NM\_000397). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113312(CYBB antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human skin tissue slide using Catalog No:113312(CYBB Antibody) at dilution of 1:200 (under 40x lens).